

Advanced  
Solutions  
in Acute  
Respiratory  
Care

# SealFlex™


Complete line of  
ribbed cushion,  
"full face" masks



SealFlex masks use the positive pressure of the ventilator to make an airtight seal, eliminating the need for excessive external pressure to secure the mask.


Just position the mask with your selected head strap with a light amount of pressure and let SealFlex do the rest.

- Double port in 3 adult sizes with pressure monitoring port
- Single port in 6 sizes from infant to large adult
- Patented design
- Standard ISO connections
- Mask can never be deflated
- Size differentiated by color for easy identification



# SealFlex™ Face Masks

The SealFlex family of products offers a combination of single and double port ribbed cushion masks.



The single and double port SealFlex offers the clinician the widest array of masks available in the market today. This highly flexible ribbed cushion design insures even the most challenging facial features can be easily sealed with minimal effort and maximum comfort for the patient. For continuous use, SealFlex promotes higher compliance for ongoing therapy.

- Soft cushion material provides greater patient comfort and tolerance
- Anatomically correct ribbed cushion shape allows for a low pressure seal, minimizing air leakage
- High dome eliminates nasal trauma
- In clinical use for over a decade

## Product Specifications

Connectors:	Double Port (Adult)
	Inlet ISO (22mm M), Outlet ISO (22mm F)
Materials:	Extended Use-TPE Cushion
Connector:	Single Port (Adult/Child)
	Inlet ISO (22mm F)
	Single Port (Infant/Neonatal)
	Inlet ISO (15mm M)
Materials:	Single Use-PVC Cushion

## Product Descriptions

Part No.:	Type:	Size:	Qty/Box:
8-209000-00	Double Port (Extended Use)	Adult-Large	10
8-209001-00	Double Port (Extended Use)	Adult-Medium	10
8-209002-00	Double Port (Extended Use)	Adult-Small	10
4-005514-00	Single Port (Extended Use)	Adult-Large	20
4-005515-00	Single Port (Extended Use)	Adult-Medium	20
4-005516-00	Single Port (Extended Use)	Adult-Small	20
4-005517-00	Single Port (Extended Use)	Child	20
4-005518-00	Single Port (Extended Use)	Infant	20
4-005519-00	Single Port (Extended Use)	Neonatal	20
8-205504-00	Single Port (Single Use)	Adult-Large	20
8-205505-00	Single Port (Single Use)	Adult-Medium	20
8-205506-00	Single Port (Single Use)	Adult-Small	20
8-205507-00	Single Port (Single Use)	Child	20
8-205508-00	Single Port (Single Use)	Infant	20
8-205509-00	Single Port (Single Use)	Neonatal	20

## Accessories

8-209102-00	Head Strap-"H"-Dry Natural Rubber	10
8-209103-00	Head Strap-"H"-Latex Free	10
4-004999-00	Head Strap-"H"-Silicone	10
8-209104-00	Head Strap-"Band"-Dry Natural Rubber	10

## Caution:

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.



Caradyne is a registered trademark.  
SealFlex is a trademark of Caradyne, Ltd.



Manufactured by:  
**Caradyne**  
Parkmore Business Centre  
Parkmore West, Galway, Ireland  
T: + 353 (0) 91 709010  
F: + 353 (0) 91 758929  
Email: info@caradyne.com  
Web: www.caradyne.com

**Caradyne**  
2055 Executive Drive  
Indianapolis, IN 46241 USA  
T: 317.246.9010  
F: 317.246.9011