

Advanced  
Solutions  
in Acute  
Respiratory  
Care

# WhisperFlow 2™

Variable  
High Flow  
Generator




The WhisperFlow 2 is the first in a series of "integrated" Continuous Positive Airway Pressure (CPAP) devices from Caradyne. WhisperFlow 2 has the same broad operating specifications as the original WhisperFlow, (i.e.) variable flow rate from 0 to 140 L/min and oxygen concentration from 28 to 100%. In addition, the WhisperFlow 2 integrates some features of a typical CPAP circuit into the generator, making set-up easier and saving money by reducing components. When used with the Criterion monitors, the caregiver has a truly integrated CPAP system with sophisticated patient monitoring. The WhisperFlow 2-40 integrates the Criterion 40 pressure monitor and the WhisperFlow 2-60 integrates the Criterion 60 pressure and oxygen monitor.

- WhisperFlow 2 generators are designed to be easily mounted to any one of our series of Criterion monitors
- Variable flow generator allows control of O<sub>2</sub> concentrations from 28 to 100% at flows from 0 to 140 L/min
- Built-in pressure relief valve at 28cm H<sub>2</sub>O
- Oxygen sensor port replaces a sensor tee piece
- Integrated pressure monitoring port
- Pole mountable for ease of transport
- Swivel elbow inlet port simplifies circuit configuration
- Two year warranty



WhisperFlow 2 and Criterion Monitor create an integrated CPAP treatment system




# WhisperFlow 2™

WhisperFlow 2 is a fully integrated CPAP solution. This generator combines traditional high performance with added features to improve ease of use.



**WhisperFlow 2**  
Variable Flow and O<sub>2</sub> Output



The WhisperFlow 2 high flow CPAP generator provides a safe, efficient solution for treating various respiratory conditions.

#### Clinical Benefits of Mask CPAP:

- Reduces patient's work of breathing
- Improves O<sub>2</sub> uptake
- Alleviates tachypnea
- Pre-emptive treatment of cardiogenic pulmonary edema
- Prevents post-operative atelectasis
- Assists in the prevention and treatment of ARDS
- Accelerates patient weaning from mechanical ventilation
- Frees up ICU ventilators in CPAP applications

#### Product Specifications

Oxygen fittings:	DISS, 6mm Hose Barb, Shrader, Mini AGA
Flow rate range:	0-140 L/min
FiO <sub>2</sub> range:	28-100%
Criterion Monitors:	See Monitor Sales Sheets

#### WhisperFlow 2 Generators

Part No:	Type:	Qty/Box:
8-100474-XX*	Europe	1
8-100478-00	North and South America	1

#### WhisperFlow 2-40

(Includes WhisperFlow 2 generator, Criterion 40 pressure monitor and pole mounting bracket)

8-100532-XX*	Europe	1
8-100533-00	North and South America	1

#### WhisperFlow 2-60

(Includes WhisperFlow 2 generator, Criterion 60 pressure and oxygen monitor and pole mounting bracket)

8-100522-XX*	Europe	1
8-100526-00	North and South America	1

#### Accessories

Part No:	Type:	Qty/Box:
8-100315-00	Pole/Rail Mounting Kit	1
8-100343-00	Caradyne Oxygen Sensor	1
8-100380-00	Oxygen Hose (84" w/DISS Fitting)	1
8-100086-00	Oxygen Hose (Schrader/Aer Liquide)	1
8-100536-00	Oxygen Hose (Schrader M/Schrader F)	1
8-100471-00	WhisperFlow 2 Filter	20
8-100264-00	Adaptor (22mm F/30mm F)	50
4-005053-00	Adaptor (22mm M/30mm M)	50
4-005054-00	Tee Piece (15mm F/22mm M), 16" Tube, Adaptor (22mm M/30mm M)	10

#### Notes:

\*Part number extensions are country specific.

#### Caution:

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.



Caradyne is a registered trademark.  
WhisperFlow is a trademark of Caradyne, Ltd.



**Caradyne**  
Parkmore Business Centre  
Parkmore West, Galway, Ireland  
T: + 353 (0) 91 709010  
F: + 353 (0) 91 758929  
Email: info@caradyne.com  
Web: www.caradyne.com

**Caradyne**  
2055 Executive Drive  
Indianapolis, IN 46241 USA  
T: 317.246.9010  
F: 317.246.9011