

Advanced  
Solutions  
in Acute  
Respiratory  
Care



# WhisperPaks™


Pre-packaged  
CPAP and  
respiratory  
treatment kits

WhisperPaks are complete compact procedure packs that contain all the necessary components to deliver therapeutic CPAP to those in need. All these packs are specifically designed to be used with the WhisperFlow high flow generator.

The WhisperPak provides over a dozen configurations to choose from. This allows you to reduce inventory requirements and training costs while establishing a CPAP standard of care.

Caradyme products that can be part of your WhisperPak are:

- Isobaric CPAP Valves
- AirSeal Face Masks
- SealFlex Face Masks
- Head Straps
- Ventilator Tubing
- Safety Valves
- Heated Wire Circuits
- CPAP Filters
- Tee Pieces



# WhisperPaks™

Part No.	Description.	Qty/Box	Part No.	Description.	Qty/Box
<b>CPAP WhisperPaks</b>			<b>Heated Circuit WhisperPaks</b>		
8-100071-00	2.5cm H <sub>2</sub> O Isobaric CPAP Valve	10	8-100377-00	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, Heated Wire Circuit	10
	2.5cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask (Adult-Large), 80" (2.0m) Tubing (22mm, Single Use)			5.0cm H <sub>2</sub> O Isobaric CPAP Valve, 20.0cm H <sub>2</sub> O Isobaric CPAP Valve, Heated Wire CPAP Circuit, Tee Piece/Safety Valve, SealFlex Double Port Mask (Adult-Medium), CPAP Filter	
8-100072-00	5.0cm H <sub>2</sub> O Isobaric CPAP Valve	10	<b>ET Tube WhisperPaks</b>		
	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask (Adult-Large), 80" (2.0m) Tubing (22mm, Single Use)		8-100356-00	ET Tube	10
8-100073-00	7.5cm H <sub>2</sub> O Isobaric CPAP Valve	10		48" (1.2m) Tubing (22mm, Single Use) (2), 16" (0.4m) Tubing (22mm, Single Use), Tee Piece/O <sub>2</sub> Analyzer (2), Tee Piece/Safety Valve, Straight Connector (22mm F/22mm F)	
	7.5cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask (Adult-Large), 80" (2.0m) Tubing (22mm, Single Use)		<b>Mask WhisperPaks</b>		
8-100074-00	10.0cm H <sub>2</sub> O Isobaric CPAP Valve	10	8-100269-00	AirSeal Double Port Mask	10
	10.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask (Adult-Large), 80" (2.0m) Tubing (22mm, Single Use)			AirSeal Double Port Mask (Adult-Large), Head Strap-Dry Natural Rubber, 80" (2.0m) Tubing (22mm, Single Use), Tee Piece/O <sub>2</sub> Analyzer, Tee Piece/Safety Valve, Straight Connector (22mm F/22mm F)	
8-100075-00	2.5cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube	10	8-100337-00	SealFlex Double Port Mask	10
	2.5cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube, 80" (2.0m) Tubing (22mm, Single Use), 0.4m Tubing (22mm, Single Use)			SealFlex Double Port Mask (Adult-Large), Head Strap-Dry Natural Rubber, 80" (2.0m) Tubing (22mm, Single Use), Tee Piece/O <sub>2</sub> Analyzer, Tee Piece/Safety Valve, Straight Connector (22mm F/22mm F)	
8-100076-00	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube	10	<b>Custom Outfitted WhisperPaks</b>		
	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube, 80" (2.0m) Tubing (22mm, Single Use), 16" (0.4m) Tubing (22mm, Single Use)		Each WhisperPak contains 1 each of all items listed unless otherwise specified.		
8-100077-00	7.5cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube	10	Custom configured WhisperPaks can be provided to meet your precise needs.		
	7.5cm H <sub>2</sub> O Isobaric CPAP Valve, Tee Piece for ET/Trach. Tube, 80" (2.0m) Tubing (22mm, Single Use), 16" (0.4m) Tubing (22mm, Single Use)				
8-100384-00	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Single Port Mask, Safety Valve	10			
	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, 20.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Single Port Mask (Adult-Medium), 80" (2.0m) Tubing (22mm, Single Use), 16" (0.4m) Tubing (22mm, Single Use), Tee Piece/Safety Valve, Tee Piece/O <sub>2</sub> Analyzer, CPAP Filter				
8-100387-00	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask, Safety Valve	10			
	5.0cm H <sub>2</sub> O Isobaric CPAP Valve, 20.0cm H <sub>2</sub> O Isobaric CPAP Valve, SealFlex Double Port Mask (Adult-Medium), 80" (2.0m) Tubing (22mm, Single Use), Tee Piece/Safety Valve, CPAP Filter				

**Caution:**

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.



0344

Caradyne is a registered trademark. WhisperPak, AirSeal and SealFlex are trademarks of Caradyne, Ltd.



Manufactured by:  
**Caradyne**  
 Parkmore Business Centre  
 Parkmore West, Galway, Ireland  
 T: + 353 (0) 91 709010  
 F: + 353 (0) 91 758929  
 Email: info@caradyne.com  
 Web: www.caradyne.com

**Caradyne**  
 2055 Executive Drive  
 Indianapolis, IN 46241 USA  
 T: 317.246.9010  
 F: 317.246.9011